


Use Instructions for 3D printed N95 Mask Clip


Purpose of clip:
to ease the placement and removal of masks and to extend the life of the mask straps



1. Cut the bottom strap of the mask



2. Place each separated strap along the clips
(adjust tightness of the mask by moving strap along the clip)



3. Clamp each strap by closing the clips




4. Put the mask on and pull the top strap over your head.



5. Grab each clip on the bottom strap and attach them together behind your neck*

TIP: if you are having trouble attaching the clips together behind your head, pinch the base of each clip with your index finger and thumb, while feeling the ball and socket of each clip with your middle finger.



*For questions or feedback, please contact Tate Park at:
tpark21@germantownfriends.org*